

## Disposable Expiratory Valve

Make your life easier. The disposable expiratory valve is the alternative to the reusable expiratory valve used for the Evita ventilator family.



It assures the functionality of intensive care ventilation and helps keep bacterial infection to a minimum.

The benefits are:

- Single-patient use
- No sterilization required
- Helps to ensure hygienic standards
- Easy handling
- Tested Dräger quality

### OPTIMIZED CLINICAL ROUTINES

In today's clinical environment, the number of multi-resistant pathogens is steadily increasing.<sup>1</sup> This type of contamination can result in time-consuming and costly consequences. Dräger's disposable expiratory valve helps to maintain hygienic standards since the valve is changed and disposed of with every patient. This is another example of how accessories may support clinical and financial benefits across all applications and how Dräger is working to "make your life easier."

### TESTED PERFORMANCE ACCORDING TO DRÄGER STANDARDS

Test procedures at Dräger help to maintain performance and safety standards needed to protect the patient and the caregiver – that is what the disposable expiratory valve does. Field tests have demonstrated that the disposable expiratory valve performs as well as conventional reusable valves and is easier to use.

### FIELD OF APPLICATION

Disposable expiratory valve for use with the Evita ventilator family (Evita 2, Evita 2 dura, Evita 4, Evita XL). For adult ventilation only.



Disposable expiratory valve

<sup>1</sup> Haley RW, Culver DH, White JW, Morgan WM, Emori TG. The nationwide nosocomial infection rate: a new need for vital statistics. Am J Epidemiol 1985;121:159-67

**TECHNICAL DATA****Material**

Housing	PP
Cover and connector	PE-HD
Sieve	Stainless steel
Screen holder	TPE
Membrane	Silicone, PE-HD, brass

**Storage**

Temperature	-20 to 60 °C
Air pressure	500 to 1060 hPa
Rel. humidity	5 to 95 %



MT-9362-2006

Disposable expiratory valve for the Evita family

**ORDER INFORMATION****Description**

Disposable expiratory valve (set of 10)

**Order no.**

8414776

**HEADQUARTERS**

Dräger Medical AG & Co. KG  
 Moisinger Allee 53–55  
 23558 Lübeck, Germany

[www.draeger.com](http://www.draeger.com)

**EUROPE, MIDDLE EAST, AFRICA,  
 LATIN AMERICA, ASIA, PACIFIC**  
 Dräger Medical AG & Co. KG  
 Moisinger Allee 53–55  
 23558 Lübeck, Germany  
 Tel +49 451 882-0  
 Fax +49 451 882-2080  
[info@draeger.com](mailto:info@draeger.com)

**CANADA**  
 Draeger Medical Canada Inc.  
 120 East Beaver Creek Road Suite 104  
 Richmond Hill Ontario L4B 4V1, Canada  
 Tel +1 905 763-3702  
 Toll-free +1 866 343 22 73  
 Fax +1 905 763-1890  
[Canada.Support@draeger.com](mailto:Canada.Support@draeger.com)

**USA**  
 Draeger Medical, Inc.  
 3135 Quarry Road  
 Telford, PA 18969, USA  
 Tel +1 215 721 5400  
 Toll-free +1 800 437 24 37  
 Fax +1 215 723 59 35  
[info.usa@draeger.com](mailto:info.usa@draeger.com)

**Manufacturer:**  
 Dräger Medical AG & Co. KG  
 23542 Lübeck, Germany  
 The quality management system at  
 Dräger Medical AG & Co. KG is  
 certified according to ISO 13485,  
 ISO 9001 and Annex II.3 of Directive  
 93/42/EEC (Medical devices).